



# DISPOSABLE ENDOSCOPY VALVES

## STERILE VALVE SETS FOR MORE EFFICIENCY AND SAFETY

The possibility of cross contamination cannot be eliminated even when preparation is done under the highest standards. The MICRO-TECH disposable valves, provide a safe and efficient solution to ensure the protection for users and patients. Time-consuming manual cleaning processes are no longer necessary. What that means for you: Simple processes, reducing time and costs, elimi-

nating the need to document reprocessing steps and reducing the risk of infection for your patients. All MICRO-TECH valves are manufactured from high quality materials and are precisely fabricated to ensure optimum functionality. The range includes valves for the most common endoscopes from Olympus, Fujifilm and Pentax.

## SPECIFIC CHARACTERISTICS

- Single use sterile packaged
- Ready for instant use
- Reduced risk of infection for users and patients
- No time consuming preparation
- More efficiency and budget security
- Compatible with endoscopes from Olympus, Fujifilm and Pentax



3-part sets (air-water-, purge and biopsy valve)

## SPECIFICATIONS

REF	Description	Package unit
<b>JET CHANNEL ADAPTER</b>		
AC03-221.F7	for Fujifilm 700 series endoscopes	Box of 100, carton 800 pieces: 8 pieces of 100
AC03-221.O	for Olympus endoscopes	Box of 100, carton 800 pieces: 8 pieces of 100
AC03-221.P	for Pentax endoscopes	Box of 100, carton 800 pieces: 8 pieces of 100
<b>ENDOSCOPE VALVE SET</b>		
AC03-212.F7	for Fujifilm 700 series with biopsy-, air-water-, suction valve	Box of 25, carton 100 pieces: 4 pieces of 25
AC03-212.O	for Olympus with biopsy-, air-water-, suction valve	Box of 25, carton 100 pieces: 4 pieces of 25
AC03-211.F7	for Fujifilm 700 series with air-water-, suction valve	Box of 25, carton 100 pieces: 4 pieces of 25
AC03-212.F	for Fujifilm 500/600 series with biopsy-, air-water-, suction valve	Box of 25, carton 100 pieces: 4 pieces of 25
AC03-213.O	for Olympus with biopsy-, air-water-, suction valve, Jet Adapter	Box of 25, carton 100 pieces: 4 pieces of 25
AC03-211.O	for Olympus with air-water-, suction valve	Box of 25, carton 100 pieces: 4 pieces of 25
<b>CLEANING VALVE</b>		
AC03-220.F7	for Fujifilm 700 series endoscopes	Box of 100, carton 400 pieces: 4 pieces of 100
AC03-220.O	for Olympus endoscopes	Box of 100, carton 400 pieces: 4 pieces of 100
<b>PROTECTIVE CAP DISTAL TIP</b>		
AC06-201	External diameter endoscope: 7-15 mm	100 pieces
AC06-202	External diameter endoscope: 2-10 mm	100 pieces