



# RETRIEVAL DEVICES

FORCEPS, BASKETS AND NETS



# FOREIGN BODY FORCEPS

## QUALITY CRAFTED FOR HIGH PRECISION CONTROL

Due to the ergonomic 3-ring handle the one-time foreign body forceps can be used with ease and comfort. Each movement is passed exactly to the forceps – even in angled body regions. The calipers themselves are made

of high quality material and are perfectly manufactured. These qualities allow you to grab firmly, guaranteeing safe foreign body retrievals. The range includes forceps for gastroenterology, pneumology, paediatrics and urology.

### SPECIFIC CHARACTERISTICS

- High-quality processing
- Covering with markings
- Ergonomic 3-ring handle
- Easy and safe operation
- Different areas of application
- 3 types: Alligator, 2-1 teeth and Griffin
- Sterile and individually wrapped
- Type Griffin also for stent-extraction and -repositioning



Type Griffin: long alligator jaw with 2:1 teeth



Type: 2-1 teeth

## SPECIFICATIONS

REF	Type	Length mm	Ø mm
<b>DISPOSABLE FOREIGN BODY FORCEPS - GASTROENTEROLOGY</b>			
NGF03-11023230	2-1 teeth	2300	2.3
NGF13-11023230	Alligator	2300	2.3
NGF33-11023230	Griffin	2300	2.3
<b>DISPOSABLE FOREIGN BODY FORCEPS - PNEUMOLOGY, PAEDIATRICS, UROLOGY</b>			
NGF32-11018180	Griffin	1800	1.8
NGF32-11018120	Griffin	1200	1.8

Packaging unit: 10 pieces

Type Griffin is ideal and our recommendation for stent-extraction and -repositioning.



# BASKET EXTRACTOR


## REMOVE FOREIGN OBJECTS SAFELY AND EFFECTIVELY

The extraction baskets from Bauer medical are not only suitable for stone extraction in the bile duct, but also for use in the upper and lower intestinal tract for the safe recovery

and removal of foreign bodies such as marbles, batteries or mouth fitting toys. The rotation of the Vario basket ensures safe and reliable recovery of the foreign body.

### SPECIFIC CHARACTERISTICS

- Fixed pre-determined breaking point
- Smooth wire with optimized gliding quality
- Ergonomic grip for improved tactile feeling
- Stone extractor basket with atraumatic tip
- Vario basket



Atraumatic tip gentle to working channel and mucosa



Highly mobile grip for precision insertion and extraction

## SPECIFICATIONS

REF	Form	Features	Ø basket mm
<b>DISPOSABLE BASKETS, PULL ROPE-LENGTH 200 CM</b>			
BK1-B10-25-23-200	Vario	Multifilament, 150°-rotation, protection head	15 / 25
BK1-B10-30-23-200	Vario	Multifilament, 150°-rotation, protection head	20 / 30
BK1-B20-20-23-200	Hexagonal	Multifilament, 4 rope, protection head	20
BK1-B20-25-23-200	Hexagonal	Multifilament, 4 rope, protection head	25
BK1-B20-30-23-200	Hexagonal	Multifilament, 4 rope, protection head	30
BK1-B20-35-23-200	Hexagonal	Multifilament, 4 rope, protection head	35
BK1-B21-35-23-200	Hexagonal	4 rope, half coiled	35
BK1-B22-20-23-200	Hexagonal	4 rope, full coiled	20
BK1-B22-25-23-200	Hexagonal	4 rope, full coiled	25
BK1-B22-30-23-200	Hexagonal	4 rope, full coiled	30
BK1-B22-35-23-200	Hexagonal	4 rope, full coiled	35
<b>DISPOSABLE GUIDE WIRE BASKETS, PULL ROPE-LENGTH 200 CM</b>			
BK1-B90-25-27-200	Combi-hexagonal	Multifilament, retractable tip	15 / 25
<b>DISPOSABLE LITHOTRIPSY BASKETS, PULL ROPE-LENGTH 400 CM</b>			
LT1-B10-25-23-400	Vario	Multifilament, 150°-rotation, retractable tip	15 / 25
LT1-B10-30-23-400	Vario	Multifilament, 150°-rotation, retractable tip	20 / 30
LT1-B20-20-23-400	Hexagonal	Multifilament, retractable tip	20
LT1-B20-25-23-400	Hexagonal	Multifilament, retractable tip	25
LT1-B20-30-23-400	Hexagonal	Multifilament, retractable tip	30
LT1-B20-35-23-400	Hexagonal	Multifilament, retractable tip	35

Packaging unit: 5 pieces



# NETS

## SAFELY REMOVE FOREIGN BODIES AND POLYPS

The MT net has been developed for the removal of polyps and foreign bodies and enables the reliable recovery of tissue samples. The net itself consists of a reinforced thread

material, which offers additional safety during recovery. Thanks to its unique dimensional stability and large opening, it facilitates work over several passes.

### SPECIFIC CHARACTERISTICS

- Large net
- High quality manufacture
- Comfortable 3-ring handle
- Easy and precise opening



Shape stable MT net

## SPECIFICATIONS

REF	Ø net mm	Snare	Working length mm	Ø working channel mm
<b>RETRIEVAL NET</b>				
NT-E-30/16-P	30	Oval	1600	2.8

Packaging unit: 10 pieces

