

**BIOPSY FORCEPS**  
**DISPOSABLE FORCEPS**



# DISPOSABLE BIOPSY FORCEPS

## RELIABLE EXAMINATION AND EVALUATION OF TISSUE SAMPLES

A reliable diagnosis requires the targeted removal of biopsy specimens. The excellent cutting edge of our biopsy forceps allows you to take diagnostically conclusive tissue samples in a safe and easy manner. The forceps are

available with diameters from 1,8 mm to 2,4 mm and with a length of 120 cm. Whether they come with or without a spike, with standard or toothed jaws – all models are characterised by precise control and easy insertion.

### SPECIFIC CHARACTERISTICS

- Large tissue samples for reliable diagnoses
- Precise control even with highly angled endoscopes
- Very comfortable to use
- Wide product variety: with or without a spike, standard or toothed
- High quality riveted scissors joints
- Working channel-friendly design



Toothed jaw



Scissors joint


## SUCCESSFUL IN PRACTICAL USE

The biopsy forceps perform in practice thanks to their ease of insertion and excellent mobility. The area to be examined can therefore be reached easily and rapidly. During tissue sampling, the forceps are outstanding for their accurate control and precise cut.

Due to their particularly sharp cutting edge, the forceps allow precise removal of biopsates. The amount of effort in this case is minimal. The biopsates acquired biopsy samples are large and therefore provide sufficient tissue for diagnostically conclusive findings.



**LARGE SPAN FOR LARGE TISSUE SAMPLES**



**COATING FOR ENDOSCOPE-FRIENDLY FEED**



## ERGONOMICAL HANDLE

With the grip of the forceps, the jaw section is perfectly under control. It offers you precise guidance even with sharply angled endoscopes.

## FORCEPS AT A GLANCE

The range includes forceps of various sizes and with different features.

### TECHBITE™



Jaw Ø 1.8 mm,  
120 cm length



Jaw Ø 2.3 mm,  
120 cm length

### TRUBITE



Jaw Ø 1.8 mm,  
120 cm length



Jaw Ø 1.8 mm,  
120 cm length

### LARGE CAPACITY



Jaw Ø 2.4 mm,  
120 cm length



Jaw Ø 2.4 mm,  
120 cm length

## SPECIFICATIONS

REF	Working channel Ø mm	Length mm	Opening width mm	Jaw Ø mm	Spike	T-Block	Covering	Jaw	Grip
<b>TECHBITE™</b>									
EBF33-11018120	2.0	1200	4.5	1.8	No	No	Grey	Smooth	1-ring
EBF33-11023120	2.8	1200	6.7	2.3	No	Yes	Grey	Smooth	1-ring
EBF33-11123120	2.8	1200	6.7	2.3	Yes	No	Grey	Smooth	1-ring
EBF43-11018120	2.0	1200	4.5	1.8	No	No	Grey	Serrated	1-ring
EBF43-11118120	2.0	1200	4.5	1.8	Yes	No	Grey	Serrated	1-ring
EBF43-11023120	2.8	1200	6.7	2.3	No	No	Grey	Serrated	1-ring
EBF43-11123120	2.8	1200	6.7	2.3	Yes	No	Grey	Serrated	1-ring
<b>TRUBITE</b>									
EBF53-11018120	2.0	1200	4.5	1.8	No	No	Grey	Smooth	1-ring
EBF53-11118120	2.0	1200	4.5	1.8	Yes	No	Grey	Smooth	1-ring
EBF13-11018120	2.0	1200	4.5	1.8	No	No	Grey	Serrated	1-ring
EBF13-11118120	2.0	1200	4.5	1.8	Yes	No	Grey	Serrated	1-ring
<b>LARGE CAPACITY</b>									
LCEBF63-10024120	2.8	1200	9.0	2.4	No	No	Grey	Serrated	1-ring
LCEBF63-10124120	2.8	1200	9.0	2.4	Yes	No	Grey	Serrated	1-ring

Packaging unit: 10 pieces

# NOTES





