



COLDSNARE

RELIABLE RESECTION OF POLYPS $\leq 9\text{MM}$

The Coldsnare from MICRO-TECH is an instrument that is suitable above all for cold resection of polyps $\leq 9\text{mm}$. This thin, braided cutting wire and its 0.23 mm wire thickness were specially developed for cold resection and makes for a highly precise, clean cut in combination


with the snare design optimized for excision of small polyps. The excised polyp is free of thermal defects and ensures that the histological assessment will provide valuable information.

SPECIFIC CHARACTERISTICS

- Compliant with guidelines for polyps $\leq 9\text{mm}$
- Special cutting wire
- Optimized snare design
- Precise, uniform cut
- High level of control
- Ergonomic grip



Precision grip with ergonomic design



SPECIFICATIONS

REF	Tube Ø mm	Length mm	Snare Ø mm	Snare type	Wire Ø mm
COLDSNARE					
CS3-11023230	2.3	2300	10	Rhombus	0.23
CS3-11523230	2.3	2300	15	Rhombus	0.23

Packaging unit: 10 pieces

CE
1639